


ILLUSTRATIE 1

3	Stable and transient transfection of plasmid DNA in eukaryotic cells.....	7
3.1	Preparation of the transfected plasmid DNA.....	9
3.1.1	<i>Workflow</i>	9
3.1.2	<i>Material.....</i>	9
3.1.3	<i>Method.....</i>	9
3.1.4	<i>Control of quality of plasmid by agarose gelectrophoresis.....</i>	13
3.1.5	<i>Spectrophotometric analysis to determine DNA concentration.....</i>	13
3.1.5.1	<i>Materials</i>	14
3.1.5.2	<i>Methods.....</i>	14
3.2	Subculturing eukaryotic cells and application of transfection methods.....	16
3.2.1	<i>Subculturing adherent cells for transfections</i>	16
3.2.1.1	<i>Material</i>	16
3.2.1.2	<i>Method</i>	17
3.2.2	<i>Scheme: cells, transfections and plasmids.....</i>	19
3.2.3	<i>Protocol for transient transfection using the calcium phosphate transfection method</i>	21
3.2.3.1	<i>Principle</i>	21
3.2.3.2	<i>Material</i>	21
3.2.3.3	<i>Method</i>	22
3.2.4	<i>Protocol for transient transfection by means of lipofectAMINE™ 2000.....</i>	23
3.2.4.1	<i>Principle</i>	23

ILLUSTRATIE 2

Powerpoint-presentatie die hoort bij syllabus (onderdeel 3.2.3)

3.2.3 Protocol for transient transfection using calcium transfection method



Phosphate buffer
Adenovirus particles
Adenovirus DNA
Add CaCl_2
DNA
Add DNA to cells
4-16 hr 37°C
Add fresh medium to cells
Remove DNA solution
7 days 37°C
Transfer medium
Infectious virus in medium

3.2.3 Protocol for transient transfection using calcium transfection method

- Ca^{2+} + phosphate \rightarrow calcium phosphate-precipitate also DNA
- on cell monolayer
- by **endocytosis**
- DNA in cytoplasm, migration to nucleus
- transient: DNA not integrated in genome \rightarrow episomal transfection
- expression analysis after 24-48 h

Multimodale input: afbeelding met kernwoorden ter ondersteuning van de tekst

Visualisatie: benadrukken kernbegrippen / structuur

Kern lesinhoud

3.2.3 Protocol for transient transfection using calcium transfection method

Protocol: see course

DNA IN COOLING BLOCK

Plate 1

HEK-cells + β -galactosidase	HEK-cells - β -galactosidase	NMepPGK β geo (+ control)
HEK-cells + β -galactosidase	HEK-cells - β -galactosidase	NMepPGK β geo (+ control)

Plate 3

coverslip HEK-cells + GFP		
coverslip HEK-cells - GFP		

Structuur bieden

3.2.3 Protocol for transient transfection using the calcium phosphate transfection method

3.2.3.1 Principle

Structuur bieden: titels en tussentitels, met nummering

Eenvoudige zinsbouw

When calcium and phosphate are mixed, in the presence of the DNA to be transfected, a calcium phosphate precipitate is formed under the right conditions, in which the DNA is coprecipitated. This precipitate is placed on the cell monolayer (60-70% confluence), which induces the contact of the precipitate with the cell membrane endocytosis. The DNA is released in the cytoplasm and migrates to the nucleus, where transcription starts. In transient transfection, the absorbed DNA is not integrated into the genome of the cell, but is broken down after a few days or lost during cell division (episomal transfection). The expression analysis is therefore performed shortly after the transfection (24 to 48 hours after the transfection).

Ca-phosphate method is used for both stable and transient transfection.

3.2.3.2 Material

Gebruik van opsommingen om
doorlopende tekst te
doorbreken

- Growing medium
- Tips 10 – 200 – 1000 μ l
- Micropipette
- 10 ml disposable pipette
- Microcentrifuge tubes 1,5 ml
- CaCl_2 /hepes
- 0,1x TE
- 2x BS / 1x Hepes

The following solutions may be created at one time for the entire group.

CaCl_2 /hepes (10 x concentrated) – 50 ml weight

- | | | | | |
|-------------------|------|--------|---------|-------|
| • CaCl_2 | 1,25 | mol/L | M 111 | ... g |
| • Hepes | 125 | mmol/L | M 238,3 | ... g |

Adjust to pH 7,05 with NaOH.

Filter sterilization through a 0,22 μ m filter.

TE – 100 ml

	mmol/l	RMM	weight g
--	--------	-----	-------------

- | | | | |
|---------------------|-----|--------|-----|
| • Tris-HCl (pH 8,0) | 1,0 | 121,14 | ... |
| • EDTA | 0,1 | 292,25 | ... |

Filter sterilization through a 0,22 μ m filter.

BS / Hepes – 50 ml

	g/l	weight g
--	-----	-------------

- | | | |
|--|------|-----|
| • Hepes (25 mmol/L) | 5,96 | ... |
| • NaCl | 16,0 | ... |
| • KCl | 0,74 | ... |
| • $\text{Na}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$ | 0,50 | ... |
| • Dextrose (D-glucose) | 2,0 | ... |

Adjust to pH 7,05 with NaOH.

Filter sterilization through a 0,22 μ m filter.

3.2.3.3 Method

Note: In order to work scientifically correctly, wells with negative control should also be refreshed.

Cell line:

Plasmid:

1-5 µg DNA = µl DNA

Belangrijke gegevens vet
gedrukt

Stapsgewijze opbouw

- Observe the cells under the microscope and estimate the % confluent cells (.... %).
- Change the medium of the wells being transfected: remove all medium and place 1,8 ml of medium per well in the 6 well plates.
- Incubate for minimum 1 hour at 37 °C.
- Create a mixture of **1-5 µg of plasmid DNA** = µl and dilute with µl TE to 100 µl in a microcentrifuge tube of 1,5 ml.
- Add **25 µl CaCl₂/hepes**.
- Add this mixture under a sterile airflow to 125 µl of BS/hepes in a 1,5 ml microcentrifuge tube (see demonstration lecturer).
- Shake this assembly by hand for **1 min** and then centrifuge briefly afterwards.
- Add the **250 µl** to one 6-well (spread over the entire well).
- Cells are incubated until testing at 37 °C in the CO₂ incubator.

Note

For the transfection that needs to be carried out twice: double the quantity and add 240 µl per well to transform.


English words

Focus op taal:
omschrijving begrip in
het Engels +
voorbeeldzin

to precipitate	Cause (a substance) to be deposited in solid form from a solution. <i>'cell proteins were then precipitated and washed in 10% trichloroacetic acid'</i>
to estimate	Roughly calculate or judge the value, number, quantity, or extent of. <i>'Experts estimate that between 70 % and 80 % of wireless networks are insecure.'</i>

ILLUSTRATIE 3

Powerpoint-presentatie die hoort bij syllabus (onderdeel 3)



Begrippen duiden door
te visualiseren

3 Stable and transient transfection of plasmid DNA in eukaryotic cells

For the work with NMf cells (Normal Murine mammary gland – fibroblastic cells) and NME cells (Normal Murine mammary gland – epithelial cells) an L1 zone is sufficient. An L2 zone is required for working with HEK cells (Human Embryonic Kidney 293-cells). Hence, the entire lab will work as a (temporary) L2 zone. The 'Safety information guide in the lab' is brought along and the safety measures are respected.

Biosafety Level 1 (BSL-1)

BSL-1 is the basic level of protection common to most research and clinical laboratories, and is appropriate for agents that are not known to cause disease in normal, healthy humans.

Biosafety Level 2 (BSL-2)

BSL-2 is appropriate for moderate-risk agents known to cause human disease of varying severity by ingestion or through percutaneous or mucous membrane exposure. Most cell culture labs should be at least BSL-2, but the exact requirements depend upon the cell line used and the type of work conducted.

(Gibco, Cell culture basics handboek, 2016).








English words

murine	Relating to or affecting mice or related rodents
mammary gland	The milk-producing gland of women or other female mammals

Focus op taal:
omschrijving begrip in
het Engels

ILLUSTRATIE 4

Powerpoint-presentatie die hoort bij syllabus

1	Cultivate and harvest bacterial cells	 	1 min 11,000 x g
2	Cell lysis		Buffer A1 250 µl Buffer A2 250 µl Buffer A3 300 µl
3	Clarification of lysate	 	5 - 10 min 11,000 x g
4	Bind DNA	 	load supernatant 1 min 11,000 x g

Schematische
voorstelling van de
syllabus

1. Cultivate and harvest bacterial cells

Use **1-5 ml** of a saturated ***E. coli* LB culture**, pellet cells in a standard benchtop microcentrifuge for **30 sec** at **11 000 x g**. Discard the supernatant and remove as much of the liquid possible.

Note

In order to obtain more yield, the following procedure is used.

Mix *E. coli* well with the LB medium. Place the culture in a 2 ml microcentrifuge tube. Centrifuge for 1 min at full force, discard the supernatant.

Repeat this again 2 times.

Use the pellet further in step 2.


2. Cell lysis

Add **250 µl Buffer A1**. Resuspend the cell pellet completely by pipetting up and down. Make sure no cell clumps remain before addition of Buffer A2!

Attention: Check Buffer A2 for precipitated SDS prior to use. If a white precipitate is visible, warm the buffer for several minutes at 30-40 °C until precipitate is dissolved completely. Cool buffer down to room temperature (20-25 °C).

Add **250 µl Buffer A2**. Mix gently by inverting the tube **6-8 times**. Do not vortex to avoid shearing of genomic DNA. Incubate at **room temperature** for up to **5 min** or until lysate appears clear.

ILLUSTRATIE 5



<http://light-microscope.net>